
Gebrauchsanleitung
Digital-Taschen-O₂-Meter CG 867
Operating Instructions
Digital O₂ meter CG 867
Mode d'emploi
Oxygénomètre digital CG 867

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Operating Instructions

Digital O₂ meter CG 867

A. Fundamental considerations for oxygen concentration measurements in water or aqueous solution

1. The principle

The digital O₂ meter and the associated oxygen electrodes function according to the Clark principle, with platinum as cathode and silver as anode in an electrolytic cell to which a polarization voltage of 800 mV is applied. Oxygen gas present in the electrolyte is reduced to OH⁻ ions at the platinum cathode. The resulting current (1 to about 100 nA) is diffusion limited and therefore proportional to the oxygen concentration in the sample solution. This current is amplified in the O₂ meter and displayed in digital form with calibration on the basis of mg O₂/liter. The background current which flows even in the absence of oxygen, is compensated during calibration.

2. The electrode

The cathode and anode are immersed in alkaline phosphate electrolyte solution inside the electrode, to give high stability readings and short polarization times. The electrolyte compartment is separated by a PTFE film from the sample solution. Only gases, but no ions, can diffuse through this film from the sample solution into the internal electrolyte. The film thickness is about 25 µm. The electrode is constructed as 4-pole plug-head device which is connected to the oxygen meter with a special cable. The built-in NTC resistor has an electrical resistance value of 100 kOhms at 25 °C.

3. Temperature behaviour and compensation

Speed of electrode reactions depends on temperature. Therefore, for a given oxygen concentration, the electrode current becomes greater at higher temperature, and conversely. This temperature effect is compensated with the NTC resistor incorporated in the electrode. For aqueous solutions, this measuring unit therefore always gives correct readings in mg O₂/liter, independent of the actual temperature. Calculatory conversion to percentage saturation values can be made with the aid of the table given in section 8.5. The values in mg O₂/liter there specified represent 100 % oxygen saturation at the respective temperatures.

4. Side reactions and interference

All substances which can diffuse through the membrane and for which -800 mV potential suffices for polarographic reduction, will be reduced in the electrode and will give rise to a corresponding

current contribution, if they are present. Interference can be caused by ions entering the electrode through porous or mechanically damaged membranes, and by diffusion of other reactive gases apart from oxygen, e.g. Cl₂, SO₂ and H₂S. These substances react in undesired manner with the electrode. Acidic or basic gases (e.g. CO₂, NH₃) change the pH value of the electrolyte solution and thus disturb the reading, particularly when measuring small oxygen concentrations. High salt concentrations in the sample solution can falsify readings too.

B. Procedure for making measurements

5. Commencing operation of the unit

5.1 Before commencing operation for the first time, open the battery compartment on the underside of the unit and insert the provided battery (1 battery, 9 V, type IEC 6 F 22). Use only leakproof batteries for replacements.

5.2 Connect the oxygen electrode (first plug-in, then screw-in).

5.3 Switch-on the unit in mode setting "mg/l" or "°C".

5.4 Battery check: The unit incorporates an automatic battery check function. A condition of low battery voltage is indicated in the liquid crystal display (LCD).

6. Polarizing the electrode

6.1 Switch to range "mg/l" or "°C". The 800 mV polarizing voltage is applied to the electrode in both settings of the switch.

6.2 Before making measurements, check that the membrane is undamaged and that there are no gas bubbles present in the electrolyte space (see section 11 "maintenance of the oxygen electrode").

6.3 Dip the electrode into saturated sodium sulfite solution (Na₂SO₃) (zero point solution). After about 15 minutes soaking in this solution, the electrode is polarized and ready for making measurements. It is not necessary to stir the solution. The end of the polarization time is indicated by stabilized reading.

7. Temperature measurement

7.1 Set the switch to position "°C".

7.2 Dip the electrode into the sample solution which is to be measured and take the temperature reading on the liquid crystal display (accuracy ± 0.5 °C in range 0 ... 40 °C).

8. Calibrating the O₂ meter for oxygen measurements

8.1 Switch to setting "mg/l".

8.2 When the electrode has been immersed for about 15 minutes in the zero point solution for polarization and is still in this solution, turn the control knob "NULLPUNKT/ZERO" to give 00.0 mg/l display reading. When recalibrating electrodes in continuous use, the polarization time is not required because the electrode is already polarized.

8.3 Rinse the electrode with water and cautiously dry it without rubbing (use cotton wool or soft filter paper). The electrode exposed in this manner to the air reaches an equilibrium corresponding to the partial pressure of oxygen and thus to saturation oxygen concentration in water at the given temperature (air calibration). Before making sensitivity setting, measure the temperature of the air or of the air-saturated water (see section 7).

Alternatively, rinse the electrode with water and then dip it into freshly air-saturated distilled water for calibration. The advection rate must be at least 10 cm/s (stir). Air-saturated water can be prepared by shaking or vigorously stirring distilled water for 10 minutes with free access of air. Ordinary tapwater can be used for calibration in applications for which an accuracy of ± 0.2 mg O₂/l suffices.

8.4 Adjust the control knob "STEILHEIT/SLOPE" to give the display reading corresponding to the "mg O₂/l" value specified in the table below (section 8.5) for the measured temperature.

8.5 Table of oxygen saturation values
Nominal values on calibration in air or air-saturated water at 1013 mbar (= 760 mmHg)

t °C	mg O ₂ /l	t °C	mg O ₂ /l	t °C	mg O ₂ /l
0	14.2				
1	13.8	11	10.7	21	8.7
2	13.4	12	10.4	22	8.5
3	13.1	13	10.2	23	8.4
4	12.7	14	10.0	24	8.3
5	12.4	15	9.8	25	8.1
6	12.1	16	9.6	26	8.0
7	11.8	17	9.4	27	7.9
8	11.5	18	9.2	28	7.8
9	11.2	19	9.0	29	7.6
10	10.9	20	8.8	30	7.5

values from "Deutsche Einheitsverfahren zur Wasseruntersuchung" (German standard methods for water analysis) System No. G2.

8.6 Values other than those in the table above hold for making measurements in salt solutions, other solvents or in solvent mixtures with or without water as one component, due to the differing solubility of oxygen in such media.

9. O₂ measurement

9.1 Switch to setting "mg/l".

9.2 Dip the calibrated electrode into the sample solution which is to be measured and take the O₂ concentration reading in mg/l on the display. The advection rate should be at least 10 cm/s (stir the sample solution).

10. Chart recorder output

10.1 The unit has a chart recorder output (0.1 mV/digit) for the O₂ measuring range, for connecting a compensation type chart recorder. The source resistance of this chart recorder output is 2 kOhms.

11. Maintenance of the oxygen electrode

11.1 The membrane may cease to function correctly after prolonged use or mechanical damage. This will falsify the readings and the electrolyte may drain out. In such cases the membrane must be replaced.

11.2 To replace the membrane head: Unscrew the measuring head of the electrode (see Fig. 1 b) containing the membrane. Take the new measuring head, with intact membrane, fill it with electrolyte solution (see Fig. 2 a) and then **slowly** screw it onto the electrode. Avoid trapping of a gas cushion inside the electrode as far as possible. The electrode is now ready for use again.

11.3 The old measuring head can be repaired. For this purpose, unscrew it for further dismantling. Remove the old O-ring (6 x 1) and the old membrane with forceps (see Fig. 2 b). Take a new membrane from the maintenance set (25 µm PTFE film), carefully mount it with the O-ring and trim at the sides with scissors. The measuring head is fully functional again after screwing-on the sleeve. The individual steps are depicted in Figs. 2 b to 2 d.

11.4 The electrode body consisting of the electrode shaft and the plug head has almost unlimited life and requires no maintenance. After prolonged use of the electrode in solutions containing high concentrations of oxygen, the silver anode may become deactivated by deposited silver phosphate. An indication of this condition is that the electrode no longer responds correctly to calibration. In this case, screw-off the measuring head as described in section 11.2 and cautiously remove the precipitated silver phosphate from the silver electrode surface with the aid of filter paper. The roughness of the glass surface with the platinum cathode (fig. 1 b) must not be changed. If necessary, it is cleaned with filtering paper. The

electrode is then ready for operation again after refilling the measuring head with new electrolyte solution (see Fig. 2 a).

12. Points to note

- 12.1 In order to ensure continued correct functioning, never open the O₂ meter CG 867.
- 12.2 Warranty is voided in the case of unauthorized tampering inside the unit or in the case of careless or deliberate damage.
- 12.3 The very high accuracy which is possible with this unit, is realized only when using the supplied electrolyte (trisodium phosphate solution) especially when replacing the membrane.

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