GEBRAUCHSANLEITUNG
ORIGINALVERSION
OPERATING MANUAL
MODE D'EMPLOI
MANUAL DE INSTRUCCIONES

HandyLab 750 HandyLab 750 EX

SI Analytics
a xylem brand

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Operating Instructions Page 35..... 61

Important notes: Before initial operation of the HandyLab 750, 750 EX, please read and observe carefully the operating instructions. For safety reasons the HandyLab 750, 750 EX may only be used for the purposes described in these present operating instructions. Please also observe the operating instructions for the units to be connected.

All specifications in this instruction manual are guidance values which are valid at the time of printing. However, for technical or commercial reasons or in the necessity to comply with the statutory stipulations of various countries, SI Analytics may perform additions to the HandyLab 750, 750 EX without changing the described properties.

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HandyLab 750, 750 EX is a portable pH meter.

Guarantee

We provide guarantee for the device described for two years from the date of purchase. This guarantee covers manufacturing faults being discovered within the mentioned period of three years. Claim under guarantee covers only the restoration of functionality, not any further claim for damages or financial loss. Improper handling/use or illegitimate opening of the device results in loss of the guarantee rights.

To ascertain the guarantee liability, please return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.

Registered trademarks

The following names are registered trademarks. For practical reasons they are shown without trademark symbol in this manual.

- MEMOSENS® ⁽¹⁾
-  MEMOSENS® ⁽²⁾
- MemoSuite® ⁽³⁾
- Sensocheck® ⁽³⁾
- Sensoface® ⁽³⁾

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Disposal

Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

Safety guidelines (HandyLab 750, 750 EX)

Caution!

Whenever it is likely that the protection has been impaired, the device shall be made inoperative and secured against unintended operation.

The protection is likely to be impaired if, for example:

- the device shows visible damage
- the device fails to perform the intended function
- prolonged storage at temperatures above 70 °C
- severe transport stresses

In this case, a professional routine test must be performed. This test should be carried out by the manufacturer.

Caution!

Take protective measures against ESD when the sensor socket is open. At measuring temperatures from -20 to +100 °C the temperature detector can be automatically identified.

The charging time of the battery may vary when a USB port with lower power output is used. Temperatures above 50 °C may result in loss of battery capacity.

Inserting the Batteries

With four AA batteries, the HandyLab has an operating time of over 1000 h. Open the battery compartment on the rear of the device. Be sure to observe the correct polarity when inserting the batteries (see markings in the battery chamber). Close the battery compartment cover and screw it handtight.

A special lithium-ion battery suited to the battery compartment is available for the HandyLab 750. The battery is recharged through the USB port.

Note: Not available for the HandyLab 750 EX (device variant for applications in hazardous locations).

A battery icon in the display indicates the battery power level:

Icon fully filled

Batteries at full capacity

Icon partially filled

Battery capacity is sufficient

Icon empty

Battery capacity not sufficient; calibration is possible, no logging

Icon blinks

Max. 10 operating hours remaining, measurement is still possible

Caution! It is absolutely necessary to replace the batteries.

Safety guidelines (HandyLab 750 EX)

Precautions for application in hazardous locations

In hazardous areas, only the **HandyLab 750 EX** may be used. The following safety guidelines must be observed:

Warning!

- Only open the battery compartment of the HandyLab 750 EX outside the hazardous location.
 - Never try to open the device. If a repair should be required, return the device to our factory.
 - Never use the USB port within the hazardous location.
-
-

Warning!

When using the HandyLab 750 EX (device variant for applications in hazardous locations) in a hazardous location, only the battery types listed below may be used. The batteries must be from the same manufacturer and of identical type and capacity. Never use new and used batteries together.

IECEx

Batteries for Application in Hazardous Locations

Batteries (4x each)	Temp. class	Ambient temperature range
Duracell MN1500	T4	$-10\text{ °C} \leq T_a \leq +40\text{ °C}$
Energizer E91	T3	$-10\text{ °C} \leq T_a \leq +50\text{ °C}$
Power One 4106	T3	$-10\text{ °C} \leq T_a \leq +50\text{ °C}$
Panasonic Pro Power LR6	T3	$-10\text{ °C} \leq T_a \leq +50\text{ °C}$

Start-Up

Connecting a Sensor

The HandyLab 750/750EX provides several connections so that many types of sensors can be used for measurement.

Note that only **one** sensor may be connected to the meter at a time.


The meter automatically recognizes a connected Memosens sensor and switches accordingly. Memosens is signaled in the display.

Caution!

Always make sure that a sensor is connected to the meter before starting measurement.

Explanation: The analog pH input of the HandyLab is an electro-meter amplifier with an extremely high-impedance. When the sensor is not in contact with the medium or not connected to the meter, electric charges on the input can generate arbitrary, stable pH or mV values which will be shown in the display.


Connections


- a Micro USB port
- b M8, 4 pins for Memosens sensors
- c Temperature probe GND
- d Temperature probe
- e pH socket (DIN 19 262)

Memosens sensors have a **cable coupling**, which allows convenient replacement of sensors while the cable remains connected to the meter. The connecting cable is connected to socket **b** (M8, 4 pins for Memosens sensors).

Switching On the Meter

When you have connected the sensor, you can switch the meter on by pressing the **on/off** or **meas** key.


When the meter is switched on with the **on/off** key, first a self test is performed and then the calibration data and settings are displayed before the meter switches to measuring mode.

When the meter is switched on with the **meas** key, it immediately switches to measuring mode.


Depending on the connected sensor and the specific measuring task, several steps for configuration and calibration must be performed as described on the following pages.

Icons

Important information about the state of the device:

Start-Up


Keypad

The keys of the membrane keypad have a noticeable pressure point.

They have the following functions:

- on/off** Switches the meter on and displays the device and calibration data (see Start-Up)
- meas** Switches the meter on / Activates measuring mode / Stops the data logger
- cal** Starts calibration
- set** Activates configuration / Confirms entries
- clock** Displays time and date, allows setting the clock using **set**
- RCL** View stored values
- STO** Holds and saves a measured value, allows setting and starting the logger by pressing **set** (page 46)
- ▲**
▼ When this icon is displayed, you can use the arrow keys for navigation.

pH Configuration

Prior to measurement, a configuration should be performed to match the connected sensor and the desired measurement performance. Furthermore, you can select the suitable calibration method.

The following table gives you an overview. Factory settings are shown in **bold print**.

Measurement

↓ **set**


"Setup" display

Select using arrow keys, confirm by pressing **set**.

Display 1	pH x.xx / pH x.xxx / mV / (°C for analog pH only)
Display 2	OFF / date + time / date / time
CAL Timer	OFF / 1 ... 99 days
CAL	AutoCal /Manual/DATA INPUT/(ISFET-Zero)/FREE CAL
CAL POINTS	1 / 2 / 3 / 1-2-3 (for AutoCal, Manual, FREE CAL)
	-01-
	(De- NIST standard 1.679/4.006/6.865/9.180
	fault)
	-02- NIST technical 1.68/4.00/7.00/10.01/12.46
	-03- SI Analytics
	technical 2.00/4.00/7.00/10.00
	-04- Knick CaliMat 2.00/4.00/7.00/9.00/12.00
	-05- Mettler Toledo 2.00/4.01/7.00/9.21
	-06- HACH 4.01/7.00/10.01/12.00
	-07- Ciba (94) 2.06/4.00/7.00/10.00
	-08- Hamilton 2.00/4.01/7.00/10.01/12.00
	-09- Reagecon 2.00/4.00/7.00/9.00/12.00
	-10- DIN 19267 1.09/4.65/6.79/9.23/12.75
	-U1- loadable via HandyLab Pilot (User)
	OFF / 0.1h / 1h / 6h / 12h
	°C / °F
	24h / 12h
	dd.mm.yy / mm.dd.yy
	NO / YES (reset to factory settings)
	Note: All data logger entries will be deleted.
▲ BUFFER SET (AutoCal, FREE CAL)	set
Auto OFF	
Temp Unit	
Time Format	
Date Format	
Default	

- ▲ This icon prompts you to select a menu item using the arrow keys -
- ▼ the selection is confirmed by pressing **set**.

Calibrating

AutoCal Calibration

(Calibration with automatic buffer recognition)

The calibration method is selected in the configuration menu. Calibration is required to adjust the sensor to the meter. It is indispensable for achieving comparable and reproducible measurement results.

Measurement

↓ **cal**

CAL
CALIMATIC

Calibration method, the number of calibration points and the buffer set have been selected in the configuration menu.

↓

CAL 1/2/3
PRESS CAL

Dip sensor in **1st**/2nd/3rd buffer solution. It does not matter which buffer solution is taken first.

↓ **cal**

pH, mV value,
temperature

Depending on the number of calibration points, the procedure described for CAL **1/2/3** is repeated.

mV value blinks until calibration is completed, then successive display of:

CAL DATA


1/2/3 CAL POINTS

ZERO POINT

SLOPE

Then the meter switches to measuring mode.

Note: To abort calibration, you can press **meas** at any time. This will be confirmed by the display message "CAL ABORTED". Exception: When you have selected "CAL POINTS 1-2-3" and the first calibration step has been completed, the calibration process cannot be stopped any more.

DATA INPUT Calibration

(Calibration by entering known sensor values)

The calibration method is selected in the configuration menu.

Measurement

↓ **cal**

CAL
DATA INPUT

↓

ZERO POINT

Use ▲▼ to select the value for the zero point.

↓ **cal**

SLOPE

Use ▲▼ to select the value for the slope.

↓ **cal**

The calibration data will be displayed successively:

Date and time


ZERO POINT

SLOPE

Then the meter switches to measuring mode.

Note: To abort calibration, you can press **meas** at any time.

Calibrating

MANUAL Calibration

(Manual calibration)

The calibration method is selected in the configuration menu.

Measurement

↓ **cal**

CAL
MANUAL

The number of calibration points has been selected in the configuration menu.

↓

CAL 1/2/3
PRESS CAL

↓ **cal**

pH display blinks
PRESS CAL

Use ▲▼ to set the temperature-corrected pH value taken from the buffer table.

↓ **cal**

mV display blinks

↓

Depending on the number of calibration points, the procedure described above for CAL 1/2/3 is repeated.

mV value blinks until calibration is completed, then successive display of:

CAL DATA


1/2/3 CAL POINTS

ZERO POINT

SLOPE

Then the meter switches to measuring mode.

Note: To abort calibration, you can press **meas** at any time. This will be confirmed by the display message "CAL ABORTED". Exception: When you have selected "CAL POINTS 1-2-3" and the first calibration step has been completed, the calibration process cannot be stopped any more.

FREE CAL Calibration

(Free selection of calibration method)

FREE CAL calibration is selected in the configuration menu.

Measurement

↓ **cal**

CAL
AutoCal blinks

Use ▲▼ to select the required calibration method (AutoCal, DATA INPUT or MANUAL).

↓ **cal**

Perform the selected calibration (see AutoCal, DATA INPUT or MANUAL calibration).

Data Logger

The Data Logger

The meter provides a data logger. **Prior to use**, it must be configured and then activated. You can choose from the following logger types:

- DIFF (signal-controlled logging of measured variable and temperature)
- INT (time-controlled logging at a fixed interval)
- DIFF+INT (combined time- and signal-controlled logging)
- SHOT (manual logging by pressing the **STO** key)


The data logger records up to 5000 entries and saves them in a circular buffer.

Already existing entries will be overwritten.

The following data are recorded: primary value, temperature, time stamp and device status.

The HandyLab Pilot software allows convenient management of the data logger. It is always the currently selected process variable which is recorded. The "STO" icon and the memory address is displayed briefly to indicate that an entry is being saved.

Display: Icons related to the data logger

Operating Modes of the Data Logger (Logger Type)

Manual logging when logger is activated (SHOT)

In this mode, a measured value is recorded when the **STO** key is pressed.

Measurement
Logger **activated**

↓ **STO**

The measured value is saved to the address of the last recorded value + 1

Manual logging when logger is deactivated

Measurement
Logger **deactivated**

↓ **STO**

Measured value is maintained
Proposed address blinks
(address of the last recorded
value + 1)


If desired: Select start address
using ▲▼.

↓ **STO**

Measured value is saved to the desired address (e.g. for overwriting an incorrect measurement).

Interval (INT)


In this mode, the measured values are cyclically recorded.

Data Logger


Difference (DIFF)

When the delta range (process variable and/or temperature) related to the last entry is exceeded, a new entry is created and the delta range is displaced upwards or downwards by the delta value. The first entry is automatically created when the data logger is started.


Difference + Interval combined (DIFF+INT)

When the delta range related to the last DIFF entry is exceeded, a new entry is created (example: entry **A**) and the delta range is displaced upwards or downwards by the delta value. As long as the measured value remains within the delta range, logging is performed at the preset interval. The first DIFF entry is automatically created when the data logger is started.

Data Logger Menu

Logger display

Select using arrow keys, confirm by pressing **set**.

Select start address and start the data logger
Deletes all entries and starts the data logger at start address 0001
Deletes all entries
Select logger type and configure: DIFF, INT, DIFF+INT, SHOT (see table below)

Overview of data logger menu (default in bold print)

Logger type	DIFF	Delta pH / mV	OFF / pH 0.01...14.00 / pH 1.00
			OFF / 1... 1000 mV / 1 mV
		Delta °C / °F	OFF / 0.1 ... 50.0 °C / 1.0 °C
			OFF / 0.1 ... 90 °F / 1.0 °F
	INT	Interval	h:mm:ss 0:00:01 ... 9:59:59 / 0:01:00
	DIFF +INT	DIFF	See logger type DIFF
	INT	See logger type INT	
SHOT	Currently selected process variable is recorded		

Data Logger

Configuring the Data Logger

Prerequisite: The data logger is stopped (press **meas**).

Measurement

↓ **STO**

Measured value is maintained

↓ **set**

Logger: CONT blinks

↓ ▼

Logger: START blinks

↓ ▼

Logger: DEL blinks

↓ ▼

Logger: SET blinks

↓ **set**

Logger: Current logger type blinks

Select desired logger type using ▲▼:
DIFF, INT, DIFF+INT or SHOT.

↓ **set**

Select the appropriate parameters using ▲▼ and confirm each selection by pressing **set**. When configuration is finished, CONT blinks. You can start the data logger by selecting START or CONT (see page 51).

Starting the Data Logger using CONT

Prerequisite: Data logger is configured. Every time the meter has been switched off, the data logger must be restarted (exception: SHOT).

Measurement

↓ **STO**

Measured value is maintained

↓ **set**

Logger: CONT blinks

↓ **set**

Address of the last recorded value + 1 blinks
(proposed start address)

If desired: Select start address using ▲▼.

↓ **set**

The measured value is saved to the selected start address (exception: SHOT).

"... FREE MEMORY" is displayed.

"LOGGER" and "active logger type" icons are displayed.

Starting the Data Logger using START

Prerequisite: Data logger is configured. All existing entries are deleted. The start address for saving the values is 0001. Every time the meter has been switched off, the data logger must be restarted (exception: SHOT).

Measurement

↓ **STO**

Measured value is maintained

↓ **set**

Logger: CONT blinks

↓ ▼

Logger: START blinks

↓ **set**

All entries will be deleted. "5000 FREE MEMORY" is displayed.

"LOGGER" and "active logger type" icons are displayed.

Data Logger

Displaying the Logger Data

Pressing the **RCL** key displays all stored values. The HandyLab Pilot software allows convenient management of the data logger.

Measurement


↓ **RCL**

The "RCL" icon and the last recorded value is displayed.


Use ▲▼ to select the desired address. Empty memory locations will also be displayed.

↓ **RCL or meas**

Return to measurement

Example:
Measured value stored at location
0026

Example:
Empty memory location 0004

Stopping the Data Logger

You can stop the data logger at any time by pressing the **meas** key.

Measurement,
logger **activated**

↓ **meas**

Data logger is stopped. "LOGGER" and "active logger type" icons are no longer displayed. It is still possible to hold a measured value by pressing **STO** and send it to any desired address.

Clearing the Data Logger

Selecting "DEL" deletes all data records.

Measurement

↓ **STO**

Measured value is maintained

↓ **set**

Logger: CONT blinks

↓ ▼

Logger: START blinks

↓ ▼


Logger: DEL blinks
PRESS SET

↓ **set**


All stored data are deleted.
"0000 DELETED" is displayed.

Error Codes and Device Messages

Error messages are indicated as "ERROR ..." on the display. Information on the sensor condition is indicated by the "Sensoface" icon (friendly, neutral, sad) possibly accompanied by an info message ("INFO ...").

Example of an error message: ERROR 8 (identical calibration media)

Example of a "Sensoface" message: INFO 1 (cal timer expired)

Sensoface provides information on the sensor condition (maintenance request). Measurement can still be performed. After a calibration, the corresponding Sensoface icon (friendly, neutral, sad) is shown together with the calibration data. Otherwise, Sensoface is only visible in measuring mode.

The most important error messages and "Sensoface" info messages are shown on the inside of the protective cover. A complete list of messages and their meanings is provided in the following tables.

SI Analytics Made in Germany

<input type="button" value="on/off"/>	Power-on and display of device / calibration data
<input type="button" value="meas"/>	Power-on / Measuring mode / Logger stop
<input type="button" value="cal"/>	Start calibration
<input type="button" value="set"/>	View / confirm device settings
<input type="button" value="clock"/>	Time / date
<input type="button" value="clock"/> + <input type="button" value="set"/>	Set time / date

HandyLab 750/ 750 EX only

<input type="button" value="RCL"/>	View stored values
<input type="button" value="STO"/>	Save measured values
<input type="button" value="STO"/> + <input type="button" value="set"/>	Set/start logger

😊	Sensor okay
😐	Calibrate the sensor soon
😞	Calibrate or replace sensor

Error messages

1 Value out of range	10 Cal media interchanged
2 ORP value out of range	11 Sensor unstable
3 Temperature error	14 Enter time and date
4 Sensor zero	21 Sensor error
5 Sensor slope	22 Sensor conflict, disconnect one sensor
8 Identical calibration media	

See manual for other error messages.

USB
(opt.)

MEMO SENS

GND

Error Codes and Device Messages

“Sensoface” Messages

The “Sensoface” icon provides information on the sensor condition:

Sensoface Meaning

Sensor is okay

Calibrate the sensor soon

Calibrate or replace the sensor

The “neutral” and “sad” Sensoface icons are accompanied by an “INFO ...” message to give a hint to the cause of deterioration.

Sensoface

Message

Cause

INFO 1

Calibration timer

INFO 3

Sensocheck

INFO 5

Zero / Slope

INFO 6

Response time

INFO 7

ISFET: Operating point (asymmetry potential)

INFO 8

ISFET: Leakage current

INFO 9

ORP offset

Error Codes and Device Messages

Error Messages

The following error messages can be shown in the display.

Message	Cause	Remedy
 blinks	Battery empty	Replace batteries
ERROR 1	pH value out of range	Check whether the measurement conditions correspond to the adjusted measuring range.
ERROR 2	ORP value out of range	
ERROR 3	Temperature value out of range	
ERROR 4	Sensor zero point too high/low	Thoroughly rinse the sensor and recalibrate. If this does not help, replace the sensor.
ERROR 5	Sensor slope too high/low	
ERROR 8	Calibration error: Identical buffers	Use a buffer solution with a different nominal value before starting the next calibration step.
ERROR 9	Calibration error: Buffer unknown	Make sure that you use the same buffer set as configured.
ERROR 10	Cal media interchanged	Repeat calibration.
ERROR 11	Measured value unstable Drift too high	Leave the sensor in the liquid until the temperature is stable. If this does not help, replace the sensor.
ERROR 14	Time and date invalid	Set time and date
ERROR 18	Configuration invalid	Restart, reset to factory settings (Setup: DEFAULT YES), configure and calibrate. If this does not help, send in the device for repair.
ERROR 19	Factory settings error	Device defective, send it in.
ERROR 21	Sensor error (Memosens)	Connect operational Memosens sensor.
ERROR 25	Buffer distance (user defined buffer table)	Reenter the buffer table.

Specifications (HandyLab 750)

Display	Sensoface	Gives you indications about the status of the sensor, analysis of zero-point/slope, setting time, frequency of calibration. Status indication (friendly, neutral, sad)
	Display	LCD STN 7-segment display with 3 lines and icons
	Status indicators	Battery power level, logger
	Notices Keypad	Hourglass [on/off], [cal], [meas], [set], [▲], [▼], [STO], [RCL], [clock]
Diagnostics functions	Sensor data (Memosens only)	Manufacturer, sensor type, serial number, operating time
	Calibration data	Calibration date, zero, slope
	Device self-test	Automatic memory test (FLASH, EEPROM, RAM)
	Device data	Device type, software version, hardware version
Data retention	Parameters, calibration data > 10 years	
Data transfer	1x Micro USB-B for data transfer to PC	
Data logger	With up to 5000 memory locations	
Calibration data logger	Recording: manual, interval- or event-controlled	
MemoLog (Memosens only)	Up to 100 Memosens calibration protocols recordable	
	- Records displayable - Directly transferable via Memosuite (USB): manufacturer, sensor type, serial no., zero-point, slope, calibration date	
Communication	USB 2.0	
	Profile	HID, installation without driver
	Use	Data exchange and configuration via HandyLab Pilot software
Nominal operating conditions	Ambient temperature	-10 ... +55 °C
	Transport/ Storage temperature	-25 ... +70 °C
Power	Relative humidity	0 ... 95 %, short-term condensing allowed
	Power supply	Batteries 4x AA alkaline or 4x NiMH (rechargeable) or 1x Li-ion-battery
	Operating time	approx 1000 h (Alkaline)
Housing	Material	PA12 GF30 + TPE
	Protection	IP66/67 with pressure compensation
	Dimensions	approx. (132 x 156 x 30) mm
	Weight	approx. 500 g
	Licensing - Compliance - Equipment safety	
EMC	EN 61326-1 (General Requirements)	
	Emitted interference	Class B (residential area)
	Immunity to interference	Industry
RoHS conformity	EN 61326-2-3 (Particular Requirements for Transmitters)	
	According to directive 2011/65/EC	

*) User-defined

- 1) Ranges depending on Memosens sensor
- 2) According to EN 60746-1, at nominal operating conditions
- 3) ± 1 count
- 4) Plus sensor error

Specifications (HandyLab 750 EX)

Memosens pH input (ISFET)	M8 socket, 4 pins, for Memosens lab cable	
	Display ranges ¹⁾	-2,000 ... +16,000 pH -2000 ... +2000 mV -50 ... +250 °C
Memosens Redox	M8 socket, 4 pins, for Memosens lab cable	
	Display ranges ¹⁾	-2000 ... +2000 mV -50 ... +250 °C
	Sensor standardization ¹⁾	ORP calibration (zero adjustment)
pH/mV (analog)	Permissible calibration range	ΔmV (offset) -700 ... +700 mV
	Connection	pH socket, DIN 19 262 (13/4 mm)
	pH range	-2 ... 16
	Decimal places *)	2 or 3
	Input resistance	1 x 10 ¹² Ω (0 ... 35 °C)
	Input current	1 x 10 ⁻¹² A (at RT, doubles every 10 K)
	Measuring cycle	approx 1s
	Measurement error ^{2,3,4)}	< 0.01 pH, TC < 0.001 pH/K
	mV range	-1300 ... +1300
	Measuring cycle	approx 1s
Measurement error ^{2,3,4)}	< 0,1 % meas. val. + 0,3 mV, TC < 0,03 mV/K	
Temperature	2 x 4 mm dia. for integrated or separate temperature detector	
	Measuring ranges	NTC 30 kΩ -20 ... +120 °C Pt 1000 -40 ... +250 °C
	Measuring cycle	approx 1s
Measurement error ^{2,3,4)}	< 0.2 K (Tamb = 23 °C); TC < 25 ppm/K	
pH-calibration		
	Operating modes^{*)}	AutoCal Calibration with automatic buffer recognition MANUAL Manual calibration with entry of individual buffer values
AutoCal buffer sets^{*)}	DATA INPUT	Data entry of zero and slope
	-01- (Default)	NIST standard 1.679/4.006/6.865/9.180
	-02-	NIST technical 1.68/4.00/7.00/10.01/12.46
	-03-	SI Analytics technical 2.00/4.00/7.00/10.00
	-04-	Knick Calimat 2.00/4.00/7.00/9.00/12.00
	-05-	Mettler Toledo 2.00/4.01/7.00/9.21
	-06-	HACH 4.01/7.00/10.01/12.00
	-07-	Ciba (94) 2.06/4.00/7.00/10.00
	-08-	Hamilton 2.00/4.01/7.00/10.01/12.00
	-09-	Reagecon 2.00/4.00/7.00/9.00/12.00
	-10-	DIN 19267 1.09/4.65/6.79/9.23/12.75
	-U1-	loadable via HandyLab Pilot (User)
Permissible calibration range	Zero point	pH 6 ... 8
	With ISFET:	-750 ... +750 mV Operating point (asymmetry)
Calibration timer^{*)}	Slope	approx. 74 ... 104 %
	Interval 1 ... 99 days, can be switched off	

Specifications (HandyLab 750 EX)

Display	Sensoface	Gives you indications about the status of the sensor, analysis of zero-point/slope, setting time, frequency of calibration. Status indication (friendly, neutral, sad)
	Display	LCD STN 7-segment display with 3 lines and icons
	Status indicators	Battery power level, logger
	Notices Keypad	Hourglass [on/off], [cal], [meas], [set], [▲], [▼], [STO], [RCL], [clock]
Diagnostics functions	Sensor data (Memosens only) Calibration data Device self-test Device data	Manufacturer, sensor type, serial number, operating time Calibration date, zero, slope Automatic memory test (FLASH, EEPROM, RAM) Device type, software version, hardware version
Data retention	Parameters, calibration data > 10 years	
Data transfer	1x Micro USB-B for data transfer to PC Never use the USB port within the hazardous location.	
Data logger	With up to 5000 memory locations	
Calibration data logger	Recording: manual, interval- or event-controlled	
MemoLog (Memosens only)	Up to 100 Memosens calibration protocols recordable - Records displayable - Directly transferable via Memosuite (USB): manufacturer, sensor type, serial no., zero-point, slope, calibration date	
Communication	USB 2.0	
	Profile	HID, installation without driver
	Use	Data exchange and configuration via HandyLab Pilot software
Nominal operating conditions	Ambient temperature	-10 ... +55 °C
	Transport/ Storage temperature	-25 ... +70 °C
Power	Relative humidity	0 ... 95 %, short-term condensing allowed
	Power supply	Batteries 4x AA alkaline or 4x NiMH (rechargeable)
Housing	Operating time	approx 1000 h (Alkaline)
	Material	PA12 GF30 + TPE
	Protection	IP66/67 with pressure compensation
	Dimensions	approx. (132 x 156 x 30) mm
	Weight	approx. 500 g
Licensing - Compliance - Equipment safety		
EMC	EN 61326-1 (General Requirements)	
	Emitted interference	Class B (residential area)
	Immunity to interference	Industry
	EN 61326-2-3 (Particular Requirements for Transmitters)	

Specifications (HandyLab 750 EX)

Explosion protection

Europe

ATEX II 1 G Ex ia IIC T4/T3 Ga

RoHS conformity

According to directive 2011/65/EC

*) User-defined

1) Ranges depending on Memosens sensor

2) According to EN 60746-1, at nominal operating conditions

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© 2015 Xylem, Inc. Version 151109
⁽²⁾ is a registered word-/trademark of Endress + Hauser Conducta GmbH & Co. KG, Gerlingen and Knick Elektronische Messgeräte GmbH & Co. KG, Berlin.
SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries.
SI Analytics GmbH