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Operating Instructions

The Microprocessor pH Meter CG 832 can be used to measure pH values, redox voltages and temperatures. No computer programming knowledge is required to operate this unit.

On power-up this pH meter switches to the same measuring function to which it was last set before switching off. The built-in nonvolatile memory ensures that all calibration data and settings for pH measurements are retained even when the mains plug is disconnected.

The values of 15 different buffer solutions together with their temperature characteristics are stored in the CG 832 for matching to the electrode function. The pH meter incorporates automatic buffer recognition for the two last used buffer solutions.

Commencing operation

1. Please make sure that the nominal mains voltage (220 V or 110 V, stated on the type specification plate) corresponds to the nominal voltage of the local mains supply.
2. Connect the unit to the local mains voltage.
3. Connect the measuring probe, e.g. a type N 1042 A combination electrode with built-in Pt 1000 resistance thermometer (see illustration, Pos. 4 and 5).
4. Switch-on the unit (press the green illuminated rocker switch on the front side). As an automatic check, all display elements light up for 3 seconds (segment test).

pH measurement with combination pH electrode N 1042 A (with built-in temperature sensor Pt 1000)

5. Switch to "pH" function with the "mode" key. The function indicators "pH" and "auto" light up (see illustrations, Pos. 2).
6. Dip the electrode into the solution which is to be measured. After allowing sufficient time for temperature equilibration, the reading is the temperature-compensated pH value.

mV measurement

7. Switch to "mV" function with the "mode" key. The function indicator "mV" lights up (see illustration, Pos. 2).
8. Dip the electrode combination (e.g. the Pt combination electrode type Pt 62 for redox voltage measurements) into the solution which is to be measured.
9. The readings in mV appear directly on the numerical display.

Temperature measurement

10. Switch to function "C" with the "mode" key. The function indicators "C" and "auto" light up (see illustration, Pos. 2).
11. Dip the Pt 1000 resistance thermometer, e.g. Type W 2130 or the combination electrode N 1042 A (with built-in resistance thermometer), into the solution whose temperature is to be measured.
12. The temperature can be read directly on the numerical display.

Note:

If the combination electrode N 1042 A with built-in resistance thermometer is not used for pH measurement with automatic temperature compensation, but instead a separate resistance thermometer (e.g. type W 2130), then connect one 4 mm plug of the W 2130 to the Pt 1000 socket and the other 4 mm plug to the reference electrode socket (see illustration, Pos. 4 and 6).

The function indicator "man" lights up when no Pt 1000 resistance thermometer is connected.

Checking the calibration of the pH meter ("calibrate")

Calibration procedure is simplest when using the combination pH electrode type N 1042 A and the same buffer solutions as used for the last previous calibration. The buffer solutions for pH = 4.01 and pH = 6.87 (at 25 °C) have been programmed into the unit at the factory for calibration, for first-time operation of the CG 832.

The order in which the buffer solutions are used can be arbitrary. It is not necessary to use first a buffer solution having a pH value close to the electrode zero point (e.g. pH = 6.87), as otherwise required in the calibration procedure for conventional pH meters.

13. Press the "mode" key for a long time (about 2 s) (see illustration, Pos. 2) to switch to "cal" mode (calibrate). The function indicators "cal", "pH" and "auto" light up.

14. After dipping the electrode into one of the two buffer solutions, the CG 832 indicates the pH value of the recognized buffer solution (e.g. pH = 4.01) after a short waiting time in flashing mode. A brief press of the "Mode" key confirms this pH value as calibration value (see Fig., Pos. 2). (For changing the at last used buffer solution, see item A.2.)

If the indicated pH value thereafter changes slightly, which may be the case due to drift of the electrode, then correction can be made by pressing the "mode" key again (see illustration, Pos. 2).

15. Rinse the electrode with distilled water and then repeat the procedure as described above using the second buffer solution. The actual temperatures and temperature differences of the buffer solutions are automatically normalized by the CG 832 to 25 °C.

16. Terminate the calibration process by pressing the "mode" key for a long time (about 2 s) again.

17. The electrode slope (in %) and then the electrode zero point, are displayed automatically. The indicator "cal" goes dark and the indicators "pH" and "auto" light up again. The pH meter with calibrated electrode is now ready for making pH measurements.

Output for chart recorder

18. The instrument has a chart recorder output (4 mm sockets) on its rear side. The source resistance of this output is $R_i = 2 \text{ k}\Omega$. This output is intended for connecting a compensation type (servo) chart recorder with floating input circuit.

pH: 59 mV/pH unit (25 °C, slope 100 %)
mV: 1 mV/mV

Error messages

19. The slope and the zero point of the connected electrode combination are calculated during the calibration procedure. Thereby all electrodes for which the found values lie within the following ranges are accepted for use:

Slope: 90 ... 105 %

Zero point: pH = 6.0 ... 7.5

In the case of values lying outside these limits, the error messages "--1" (slope) and "--2" (zero point) are displayed. The instrument waits until the respective error message is confirmed by pressing the "Mode" key. Then it displays the out of range value in flashing mode. Finally "--" is displayed for 3 seconds. If a displacement key ("Δ" or "▽") (see illustration, Pos. 3) is pressed whilst this display is present, then the instrument will take over the calibration values even when they lie outside the specified limits.

If no key is pressed whilst the display "--" is present, then the executed calibration is rejected and the results are not stored in memory. The calibration procedure must then be repeated with fresh buffer solutions and cleaned or new electrode.

Annex: Further instructions for working with the CG 832

A1 Manual temperature compensation

When using an electrode without built-in temperature sensor (e.g. a combination pH electrode N 62), the temperature value may be set manually in the "cal" function or in function temperature measurement "C".

In this case the calibration "cal" commences with display of the manually set temperature. The displayed temperature value can be changed with the keys "Δ" and "▽" to the measuring temperature (see illustration, Pos. 3).

Thereafter briefly press the "mode" key to switch over to "pH" (see illustration, Pos. 2).

A2 Calibration with other buffer solutions

If it is the intention to use buffer solutions which are not the same as those used for the last previous calibration, then the CG 832 will not indicate a correct buffer value in the flashing display after calling the calibration routine. In this case the correct buffer can be set with the two displacement keys "Δ" and "▽" (see illustration, Pos. 3).

The buffer values which have been programmed in memory are: pH = 1.09; 1.68; 2.00; 3.56; 3.78; 4.00; 4.01; 4.65; 6.87; 7.00; 7.41; 9.18; 10.00; 10.01; 12.45 (all values for 25 °C). All temperature coefficients of these buffer solutions for the range 0 °C to 100 °C have been stored too in the memory of microprocessor.

A3 Single point calibration

For simplicity, a calibration can be made with only one buffer solution, as described in item 14. Then quit the calibration routine by pressing the "mode" key for a long time again (see illustration, Pos. 2). This single point calibration makes only a zero point correction. It does not correct for the electrode slope.

A4 Reset

Strong external interference (e.g. a lightning stroke) may disturb program execution in the unit. In this case a reset is required. For this purpose, switch the unit off, hold one of the displacement keys "Δ" or "▽" depressed and at the same time switch the unit on again (see illustration, Pos. 3). The display test then executes automatically.

The unit must be recalibrated for pH measurement after a reset.

A5 Range overshoot

If a reading is out of range (e.g. greater than 1999 mV), then the overflow condition is shown as a "--" followed by a decimal point in the display.

Safety hints

For functional and safety reasons, never open the unit. Warranty becomes void in the case of unauthorized tampering inside the unit or in the case of careless or deliberate damage.

Please also observe the instructions for electrodes.

This apparatus conforms to protection class 1.

It has been built and tested according to DIN 57 411 Part 1/WDE 0411 Part 1, protection measures for electronic measuring equipment and has left the factory in safe and proper functioning condition. In order to preserve this proper condition and to ensure safe operation, the user must obey the instructions and observe the warnings contained in this operating instructions manual.

Before switching-on, make sure that the operating voltage setting on the apparatus corresponds to the nominal value of the mains voltage.

Insert the mains plug only into a power socket which has a safety grounding contact. The protection ground line must not be interrupted by an extension cable which does not have a corresponding ground conductor.

Any interruption of the protection grounding line inside or outside the apparatus, or disconnection of the grounding conductor, can produce dangerous states of the apparatus. Deliberate disconnection or interruption of the grounding line is forbidden.

Make sure that only fuses of the specified type and with the specified nominal current rating are used as replacements. The use of repaired fuses or short circuiting of the fuseholder is forbidden.

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Schott Geräte GmbH